

Professional Oxidation Reduction Potential Electrode

ORP ELECTRODE

Model : ORP-14

ISO-9001, CE, IEC1010



Features :

- * Professional ORP electrode with high accuracy and extreme reliability.
- * Epoxy body, BNC connector and 3 meters cable.
- * 12.5 mm Dia, x 150 mm.
- * To cooperate ORP meter (ORP-203, PH-207, PH-208, YK-2001PH), application for water conditioning, aquariums, beverage, fish hatcheries, food processing, photography, laboratory paper industry, quality control.



LUTRON ELECTRONIC

The Art of Measurement